

Effortless and Streamlined Workflow for Absolute Quantitation of Therapeutic Monoclonal Antibodies using Promise Proteomics mAbXmise kits and TSQ Altis Plus Mass Spectrometer

Yvonne E. Song¹, Dorothee Lebert², Jean-Valéry Guillaubez², Stephanie Samra¹, Edward Goucher¹, Bradley Hart¹
Thermo Fisher Scientific, San Jose, CA, United States¹, Promise Proteomics, Grenoble, France²

Abstract

Purpose: To simultaneously quantify the signature peptides of multiple therapeutic monoclonal antibodies in human serum in a streamlined and confident way using Promise Proteomics kits and the TSQ Altis™ Plus mass spectrometer.

Methods: Samples were purified and trypsin digested according to the instructions provided by Promise Proteomics with minor modifications (refer to instructions manual – <https://customer.mabxmise.com>). LC-MS/MS analysis was performed using Thermo Scientific™ Vanquish™ Flex UHPLC system coupled to TSQ Altis™ Plus mass spectrometer.

Results: Successfully completed process replicate analysis and LC-MS/MS analytical performance evaluation and passed acceptance criteria according to the collaboration study plan.

Introduction

Laboratories continuously seek improved productivity and efficacy for clinical testing, ultimately impacting the turn-around time and sample throughput. Various efforts have been made including rapid and high-throughput testing, automation, and high-end instrumentation. In particular, mass spectrometry (MS) has gained significant popularity in clinical laboratories for TDM of mAbs due to its great versatility to measure such complex biological proteins qualitatively and quantitatively. Here we present the streamlined workflow for the absolute quantitation of therapeutic mAbs using the Promise Proteomics mAbXmise kits and the Thermo Scientific™ TSQ Altis™ Plus mass spectrometer.

Materials and methods

Sample Preparation

The workflow is described in Figure 1 including the plate scheme and the list of kits is below.


- ITDM1: inflammation-TDM kit 1 for quantitation of 2 mAbs including adalimumab and infliximab (mAbXmise® ITDM1, Promise Proteomics, France)

- OTDM1: oncology-TDM kit 1 for quantitation of 7 mAbs including bevacizumab, cetuximab, ipilimumab, nivolumab, pembrolizumab, rituximab, and trastuzumab (mAbXmise® OTDM1, Promise Proteomics, France)

Test Method(s)

LC-MS/MS analysis was performed using Thermo Scientific™ Vanquish™ Flex UHPLC system coupled to TSQ Altis™ Plus mass spectrometer. LC and MS conditions are listed in Table 1. Transitions can be found in Tech Note (TN001753). After sample preparation, each well was injected once for process replicate analysis. The same calibrator and QCs were pooled and redistributed for analysis of LC-MS/MS analytical performance in triplicate.

Figure 1. Workflow and plate scheme



Data Analysis

Data acquisition, processing, and reporting were performed using Thermo Scientific™ TraceFinder™ 5.1 software and the evaluation criteria for LC-MS analytical performance are listed in Table 2.

Table 1. LC and MS conditions

LC gradient			
Time (min)	% A	% B	Curve
0.0	95	5	5
0.5	95	5	5
1.0	80	20	5
4.8	60	40	5
5.5	50	50	5
6.1	10	90	5
7.5	10	90	5
7.6	95	5	5
10	95	5	5

Separation conditions	
Mobile phase A	0.1 % formic acid in water
Mobile phase B	0.1 % formic acid in 10: 10: 80 water: isopropanol: acetonitrile (v/v/v)
Flow rate	0.25 mL/min
Column temperature	40 °C (Still Air)
Injection volume	10 µL (5 µL for full-process replate analysis)

Global parameters	
Source type	Heated electrospray ionization (H-ESI)
Polarity	Positive
Spray voltage (V)	3500
Sheath gas (Arb)	40
Aux gas (Arb)	7
Sweep gas (Arb)	1
Ion transfer tube temp (°C)	325
Vaporizer temp (°C)	275
Divert valve A	0.0 min: position 1-6 (waste) 0.6 min: position 1-2 (MS) 6.2 min: position 1-6 (waste)
Probe position (x-y-z)	Center - 1.5 - LM

SRM scan parameters	
Cycle time (sec)	0.35
Q1 resolution (FWHM)	0.7
Q3 resolution (FWHM)	0.7
CID gas (mTorr)	1.5
Source fragmentation (V)	0
Chromatographic peak width (sec)	6
RF lens (V)	60
Dwell time priority	3 (normal)

Table 2. Evaluation criteria for LC-MS/MS analytical performance

Analytical characteristics	Acceptance criteria
QC accuracy	Mean concentration from 85 to 115 % for QC1 and QC2 samples
QC precision	CV ≤ 15 % for QC1 and QC2 samples
Lower limit of quantification (LLOQ)	<ul style="list-style-type: none"> The analyte mean response at the LLOQ (CAL1) > 5 times the analyte response of the zero calibrator (CAL0 = blank sample). LLOQ should be 2 µg/mL Accuracy: mean concentration from 80 to 120 % Precision: CV ≤ 20 %
Linear range	Linear fit: R ² ≥ 0.99 on the expected linear range 2 ~ 100 µg/mL

Results

The highly comparable results were observed across all six replicates for both OTDM1 and ITDM1 kits with great linearity with R² > 0.994 and % RSD < 15, indicating the reproducibility of the full plate sample preparation. Then, the six replicates of each calibrator and QC from each kit were pooled together, distributed to another plate following the same plate scheme, and analyzed in triplicates, generating 18 data points for each calibrator and QC sample. We successfully evaluated the LC-MS/MS analytical performance, showing excellent linearity with R² > 0.997 and % RSD < 10 (Figure 3 and Tables 3 and 4). Additionally, Figure 4 shows extremely reproducible retention times of all SIL-mAb peptides with RT difference ± 0.03 minutes during the evaluation.

Figure 2. Representative XICs of OTDM1 and ITDM1 peptides




Figure 3. Calibration curves for OTDM1 (blue) and ITDM1 (green) peptides (#data points=18)





Table 3. Results of LC-MS/MS analytical performance of ITDM1 peptides

Analytical characteristics	QC accuracy		QC precision		LLOQ			Linear range	
	QC1	QC2	QC1	QC2	CAL0 response (5 x CAL0)	CAL1 response	Accuracy		
Adalimumab APYTFGGGQTK	107.0%	95.2%	3.1%	1.9%	0.050 (0.250)	0.267	99.0%	3.2%	0.9984
Infliximab SINSATHYAESVK	99.8%	100.4%	8.7%	3.7%	0.002 (0.010)	0.193	111.6%	8.4%	0.9970
Infliximab SAVYLQMTDLR	103.0%	98.7%	2.6%	1.9%	0.016 (0.080)	0.219	100.2%	2.1%	0.9986
Infliximab DILLTQSPAILSVSPGER	96.4%	96.2%	2.8%	2.6%	0.024 (0.120)	0.223	105.4%	6.0%	0.9987

Table 4. Results of LC-MS/MS analytical performance of OTDM1 peptides

Analytical characteristics	QC accuracy		QC precision		LLOQ			Linear range	
	QC1	QC2	QC1	QC2	CAL0 response (5 x CAL0)	CAL1 response	Accuracy		
Bevacizumab FTFLDTSK	101.1%	102.2%	1.9%	1.6%	0.009 (0.043)	0.138	98.3%	3.8%	0.9996
Cetuximab YAESISGIPSR	97.3%	96.6%	1.5%	1.8%	0.003 (0.015)	0.243	98.3%	3.6%	0.9996
Ipilimumab GLEWTFYSDGNKK	95.0%	101.2%	6.1%	2.9%	0.009 (0.045)	0.136	106.3%	7.4%	0.9982
Nivolumab ASGITFSNSGMHWVR	95.2%	94.0%	3.3%	3.2%	0.040 (0.200)	0.297	102.6%	6.2%	0.9983
Pembrolizumab DLPLTFGGGTK	100.3%	100.0%	2.0%	1.4%	0.014 (0.070)	0.294	96.0%	1.8%	0.9996
Rituximab FSGSGGTSYSLTISR	101.9%	105.3%	6.4%	5.9%	0.002 (0.010)	0.285	100.8%	10.6%	0.9971
Trastuzumab DTYIHVWR	101.2%	101.3%	3.1%	2.3%	0.010 (0.050)	0.116	99.1%	3.9%	0.9994
Trastuzumab FTISADTSK	100.2%	99.4%	1.8%	0.8%	0.028 (0.140)	0.177	102.6%	4.4%	0.9996

Figure 4. Variation in observed retention times of all SIL-mAb peptides during the evaluation of LC-MS analytical performance



Trademarks/licensing

© 2023 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

mAbXmise monoclonal antibodies quantification kits OTDM1 and ITDM1 are In Vitro Diagnostic Medical Devices for laboratory professional use and CE-IVD labeled for Europe. Assay results are intended to be used by healthcare professionals. The kits are designed to perform absolute quantification by LC-MS (Liquid Chromatography – Mass spectrometry) of specific therapeutic monoclonal antibodies (mAbs) in a patient sample.

With respect to Materials required but not provided, the Instructions for Use state under Equipment: Triple quadrupole MS instrument or high-resolution MS instrument coupled to LC system.

PROMISE Proteomics products are distributed globally, so uses, applications, and availability of product in each country depend on local regulatory registration status. Products not registered for In Vitro Diagnostic Use are for Research Use Only. mAbXmise is a trademark of Promise Proteomics.